

盟基生技優選， 您的實驗保證！

- ✓ Targetmol 來自於美國波士頓，提供高品質研究化學品，包含 CCK-8 Kit、Protease & Phosphate inhibitor 等。
- ✓ 提供超過20,000多個小分子抑制劑及120多個化合物庫。
- ✓ 近上萬篇的學術文獻發表使用
- ✓ 每批產品都以嚴格的HNMR和HPLC技術檢測，以保障產品高純度&高質量。



Cell Counting Kit-8 (CCK-8)

非選不可的理由：

- ✓ 操作簡單，一瓶就搞定！
- ✓ 對細胞無毒性
- ✓ 高穩定度！4°C避光2年，可長時間保存於-20 °C
- ✓ 靈敏度高於其他方法，如MTT, XTT or WST-1

快速偵測：

1. 細胞存活率
2. 細胞增生
3. 細胞毒性



特性	MTT	XTT	WST-1	CCK-8
溶解度	-	+	+	+
形式	粉末	2瓶溶液	1瓶溶液	1瓶溶液
製劑	使用前溶解	使用前混合	直接使用	直接使用
靈敏度	+	++	++	+++
檢測速度	+	++	++	+++
波長	560~600nm	420~480nm	420~480nm	430~490nm
毒性	+	-	-	-
穩定度	+	-	+	++
96-well plate 兼容性	+	++	++	++
便利性	+	++	++	+++

貨號

C0005

包裝

1 mL (100 tests)

買四送一



Code: 2022Q101-02

Omics Bio
盟基生物科技

Web: www.omicsbio.com.tw
Mail: info@omicsbio.com
免付費電話: 0800-222-082

台北總公司
T: 02-8698-2268
F: 02-8698-2228


高雄分公司
T: 07-223-6080
F: 07-223-8205



保護蛋白質大作戰!

給您 **最直覺** 的選擇!

- ✓ 溶液形式, 易稀釋簡單使用 (它牌錠劑含賦形劑, 有難溶風險)
- ✓ 清楚且透明的配方材料, 可避免實驗的不確定性



買七送三

■ Protease Inhibitor Cocktail 趕走所有搗蛋壞酵素!

Ingredient	Conc. (100x)	Target	Type
AEBSF.HCl	104 mM	serine proteases	Irreversible
Aprotinin	80 μ M	serine proteases	Reversible
Bestatin	5 mM	aminopeptidase B; LAP	Reversible
E-64	1.5 mM	cysteine proteases	Irreversible
Leupeptin	2 mM	cysteine proteases; serine proteases	Reversible
Pepstatin A	1.5 mM	aspartic proteinase	Reversible

貨號	品名	包裝
C0001	Protease Inhibitor Cocktail (100X DMSO)	1 mL (DMSO)

■ Phosphatase Inhibitor Cocktail 強力盾牌給蛋白滿滿的保護力!!

Cocktail	Ingredient	Conc. (100x)	Target	Type
I (ddH ₂ O) 1mL	Sodium Fluoride	100 mM	Acid Phosphatase	Reversible
	Sodium Orthovanadate	100 mM	Alkaline phosphatase, PTPs, ATPases	Reversible
	Sodium Molybdate	115 mM	Acid and phosphoprotein Phosphatase	Irreversible
	Sodium Tartrate	400mM	Acid Phosphatase	Reversible
	Imidazole	200 mM	Alkaline phosphatase	Reversible
II (DMSO) 1mL	(-)-p-Bromotetramisole oxalate	2.5 mM	Alkaline phosphatases	Irreversible
	Cantharidin	500 μ M	Ser/Thr phosphatases	Reversible
	Microcystin LR	500 nM	Acid and PP1 and PP2A	Reversible

I+II 給您完整的保護力!

貨號	品名	包裝
C0002	Phosphatase Inhibitor Cocktail I (100 \times ddH ₂ O)	1 mL (ddH ₂ O)
C0003	Phosphatase Inhibitor Cocktail II (100 \times DMSO)	1 mL (DMSO)

兩品項
一起買
優惠 **85折!**