

Elabscience TOP100 ELISA kit

Target	Catalog	Product Name	Method	Reactivity	Detection Range	Sensitivity	定價	買一送一
8-epi-PGF2α	E-EL-0041	8-epi-PGF2α(8-Epi-Prostaglandin F2 Alpha) ELISA Kit	Competitive	Universal	15.63~1000pg/mL	9.38pg/mL	19,800	9,900
8-OHdG	E-EL-0028	8-OHdG(8-Hydroxydeoxyguanosine) ELISA Kit	Competitive	Universal	1.56~100ng/mL	0.94ng/mL	19,800	9,900
ACH	E-EL-0081	ACH(Acetylcholine) ELISA Kit	Competitive	Universal	15.63~1000pg/mL	9.38pg/mL	19,800	9,900
ACTH	E-EL-R0048	Rat ACTH(Adrenocorticotrophic Hormone) ELISA Kit	Competitive	Rat	15.63~1000pg/mL	9.38pg/mL	19,800	9,900
ADMA	E-EL-0042	ADMA(Asymmetrical Dimethylarginine) ELISA Kit	Competitive	Universal	15.63~1000ng/mL	9.38ng/mL	19,800	9,900
ADP/Acrp30	E-EL-H6122	Human ADP/Acrp30(Adiponectin) ELISA Kit	Sandwich	Human	0.16~10ng/mL	0.09ng/mL	19,800	9,900
Ang-II	E-EL-H0326	Human Ang-II(Angiotensin II) ELISA Kit	Competitive	Human	31.25~2000pg/mL	18.75pg/mL	19,800	9,900
APLN	E-EL-H0456	Human APLN(Apelin) ELISA Kit	Competitive	Human	62.5~4000pg/mL	37.5pg/mL	19,800	9,900
ApoB48	E-EL-H6084	Human ApoB48(Apolipoprotein B48) ELISA Kit	Sandwich	Human	1.56~100ng/mL	0.94ng/mL	24,360	12,180
Aβ1-40	E-EL-H0542	Human Aβ1-40(Amyloid Beta 1-40) ELISA Kit	Sandwich	Human	15.63~1000pg/mL	9.38pg/mL	19,800	9,900
BDNF	E-EL-H0010	Human BDNF(Brain Derived Neurotrophic Factor) ELISA Kit	Sandwich	Human	31.25~2000pg/mL	18.75pg/mL	19,800	9,900
C3a	E-EL-M0337	Mouse C3a(Complement Component 3a) ELISA Kit	Sandwich	Mouse	0.156~10ng/mL	0.1ng/mL	24,360	12,180
C5a	E-EL-MK1612	Monkey C5a(Complement Component 5a) ELISA Kit	Sandwich	Monkey	0.16~10ng/mL	0.09ng/mL	27,320	13,660
CPP	E-EL-H0851	Human CPP(Copeptin) ELISA Kit	Competitive	Human	31.25~2000pg/mL	18.75pg/mL	24,360	12,180
CRP	E-EL-H0043	Human CRP(C-Reactive Protein) ELISA Kit	Sandwich	Human	0.39~25ng/mL	0.23ng/mL	19,800	9,900
CTXI	E-EL-H0835	Human CTXI(Cross Linked C-telopeptide of Type I Collagen) ELISA Kit	Sandwich	Human	0.16~10ng/mL	0.1ng/mL	24,360	12,180
CTXI	E-EL-M3023	Mouse CTX I (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit	Sandwich	Mouse	31.25~2000pg/mL	18.75pg/mL	24,360	12,180
CTXI	E-EL-R1456	Rat CTXI(Cross Linked C-telopeptide of Type I Collagen) ELISA Kit	Sandwich	Rat	0.16~10ng/mL	0.1ng/mL	24,360	12,180
CTX-II	E-EL-H0837	Human CTX-II(Cross Linked C-telopeptide of Type II Collagen) ELISA Kit	Sandwich	Human	0.16~10ng/mL	0.1ng/mL	19,800	9,900
DA	E-EL-0046	DA(Dopamine) ELISA Kit	Competitive	Universal	31.25~2000pg/mL	18.75pg/mL	19,800	9,900
E2	E-EL-0150	Human/Monkey/Mouse E2(Estradiol)	Competitive	Human/Monkey/Mouse	9.37~600pg/mL	5.63pg/mL	19,800	9,900
E2	E-EL-0152	Rat/Porcine E2(Estradiol) ELISA Kit	Competitive	Rat/Porcine	1.56~100pg/mL	1.00pg/mL	19,800	9,900
EPI	E-EL-0045	EPI(Epinephrine/Adrenaline) ELISA Kit	Competitive	Universal	31.25~2000pg/mL	18.75pg/mL	19,800	9,900
ET-1	E-EL-H0064	Human ET-1(Endothelin 1) ELISA Kit	Sandwich	Human	1.25~80pg/mL	0.75pg/mL	19,800	9,900
F1+2	E-EL-H1793	Human F1+2(Prothrombin Fragment 1+2) ELISA Kit	Sandwich	Human	46.88~3000pg/mL	28.13pg/mL	24,360	12,180
FGF21	E-EL-H0074	Human FGF21(Fibroblast Growth Factor 21) ELISA Kit	Sandwich	Human	31.25~2000pg/mL	18.75pg/mL	24,360	12,180
FGF23	E-EL-H1116	Human FGF23(Fibroblast Growth Factor 23) ELISA Kit	Sandwich	Human	15.63~1000pg/mL	9.38pg/mL	24,360	12,180
FNDC5	E-EL-H2254	Human FNDC5(Fibronectin type III domain-containing protein 5) ELISA Kit	Sandwich	Human	0.16~10ng/mL	0.1ng/mL	24,360	12,180
FSH	E-EL-R0391	Rat FSH(Follicle Stimulating Hormone) ELISA Kit	Sandwich	Rat	3.13~200ng/mL	1.88ng/mL	19,800	9,900
GAL3	E-EL-H1470	Human GAL3(Galectin 3) ELISA Kit	Sandwich	Human	0.16~10ng/mL	0.1ng/mL	24,360	12,180
GFAP	E-EL-H6093	Human GFAP(Glial Fibrillary Acidic Protein) ELISA Kit	Sandwich	Human	15.63~1000pg/mL	9.38pg/mL	24,360	12,180
GHRL	E-EL-H1919	Human GHRL(Ghrelin) ELISA Kit	Competitive	Human	0.16~10ng/mL	0.1ng/mL	24,360	12,180
GLP-1	E-EL-H6025	Human GLP-1(Glucagon Like Peptide 1) ELISA Kit	Sandwich	Human	1.56~100pg/mL	0.94pg/mL	19,800	9,900
GPX1	E-EL-H5410	Human GPX1(Glutathione Peroxidase 1) ELISA Kit	Sandwich	Human	31.25~2000pg/mL	18.75pg/mL	24,360	12,180
GSH	E-EL-0026	GSH(Glutathione) ELISA Kit	Competitive	Universal	1.56~100μg/mL	0.94μg/mL	19,800	9,900
HBsAg	E-EL-H6080	Human HBsAg(Hepatitis B surface antigen) ELISA Kit	Sandwich	Human	0.31~20ng/mL	0.19ng/mL	24,360	12,180
Hepc	E-EL-H6013	Human Hepc(Hepcidin) ELISA Kit	Sandwich	Human	62.5~4000pg/mL	37.5pg/mL	19,800	9,900
HIF-1α	E-EL-H6066	Human HIF-1α(Hypoxia Inducible Factor 1 Alpha) ELISA Kit	Sandwich	Human	62.5~4000pg/mL	37.5pg/mL	19,800	9,900
HIS	E-EL-0032	HIS(Histamine) ELISA Kit	Competitive	Universal	1.56~100ng/mL	0.94ng/mL	19,800	9,900
HMGB-1	E-EL-H1554	Human HMGB-1(High Mobility Group Protein B1) ELISA Kit	Sandwich	Human	31.25~2000pg/mL	18.75pg/mL	24,360	12,180
HMGB-1	E-EL-M0676	Mouse HMGB-1(High Mobility Group Protein B1) ELISA Kit	Sandwich	Mouse	15.625~1000pg/mL	9.38pg/mL	24,360	12,180
HS	E-EL-H2364	Human HS(Heparan Sulfate) ELISA Kit	Sandwich	Human	0.31~20ng/mL	0.19ng/mL	24,360	12,180
hs-CRP	E-EL-H5134	Human hs-CRP(high-sensitivity C-Reactive Protein) ELISA Kit	Sandwich	Human	15.63~1000pg/mL	9.38pg/mL	19,800	9,900
IDO	E-EL-H2162	Human IDO(Indoleamine-2,3-Dioxygenase) ELISA Kit	Sandwich	Human	0.31~20ng/mL	0.19ng/mL	24,360	12,180

Elabscience TOP100 ELISA kit

Target	Catalog	Product Name	Method	Reactivity	Detection Range	Sensitivity	定價	買一送一
IFN-γ	E-EL-H0108	Human IFN-γ(Interferon Gamma) ELISA Kit	Sandwich	Human	15.63~1000pg/mL	9.38pg/mL	19,800	9,900
IGF-1	E-EL-H0086	Human IGF-1(Insulin-like Growth Factor 1) ELISA Kit	Sandwich	Human	1.56~100ng/mL	0.94ng/mL	24,360	12,180
IL-10	E-EL-H0103	Human IL-10(Interleukin 10) ELISA Kit	Sandwich	Human	7.81~500pg/mL	4.69pg/mL	19,800	9,900
IL-10	E-EL-R0016	Rat IL-10(Interleukin 10) ELISA Kit	Sandwich	Rat	31.25~2000pg/mL	18.75pg/mL	19,800	9,900
IL-17	E-EL-H0105	Human IL-17(Interleukin 17) ELISA Kit	Sandwich	Human	31.25~2000pg/mL	18.75pg/mL	19,800	9,900
IL-1β	E-EL-H0149	Human IL-1β(Interleukin 1 Beta) ELISA Kit	Sandwich	Human	7.81~500pg/mL	4.69pg/mL	19,800	9,900
IL-1β	E-EL-M0037	Mouse IL-1β(Interleukin 1 Beta) ELISA Kit	Sandwich	Mouse	7.81~500pg/mL	4.69pg/mL	19,800	9,900
IL-1β	E-EL-R0012	Rat IL-1β(Interleukin 1 Beta) ELISA Kit	Sandwich	Rat	31.25~2000pg/mL	18.75pg/mL	19,800	9,900
IL-23	E-EL-H0107	Human IL-23(Interleukin 23) ELISA Kit	Sandwich	Human	39.06~2500pg/mL	23.44pg/mL	24,360	12,180
IL-35	E-EL-H2443	Human IL-35(Interleukin 35) ELISA Kit	Sandwich	Human	15.63~1000pg/mL	9.38pg/mL	24,360	12,180
IL-6	E-EL-H0102	Human IL-6(Interleukin 6) ELISA Kit	Sandwich	Human	7.81~500pg/mL	4.69pg/mL	19,800	9,900
IL-6	E-EL-M0044	Mouse IL-6(Interleukin 6) ELISA Kit	Sandwich	Mouse	31.25~2000pg/mL	18.75pg/mL	19,800	9,900
IL-6	E-EL-R0015	Rat IL-6(Interleukin 6) ELISA Kit	Sandwich	Rat	12.5~800pg/mL	7.5pg/mL	19,800	9,900
IL-8	E-EL-H6008	Human IL-8(Interleukin 8) ELISA Kit	Sandwich	Human	7.81~500pg/mL	4.69pg/mL	19,800	9,900
INS	E-EL-R3034	Rat INS(Insulin) ELISA Kit	Sandwich	Rat	6.25~400pg/mL	3.75pg/mL	24,360	12,180
IP3	E-EL-0059	IP3(Inositol Triphosphate) ELISA Kit	Competitive	Universal	15.63~1000pg/mL	9.38pg/mL	19,800	9,900
KL	E-EL-H5451	Human KL(Klotho) ELISA Kit	Sandwich	Human	0.31~20ng/mL	0.19ng/mL	24,360	12,180
LEP	E-EL-H6017	Human LEP(Leptin) ELISA Kit	Sandwich	Human	15.63~1000pg/mL	9.38pg/mL	19,800	9,900
LEP	E-EL-R0582	Rat LEP(Leptin) ELISA Kit	Sandwich	Rat	0.16~10ng/mL	0.1ng/mL	19,800	9,900
LH	E-EL-M3053	Mouse LH(Luteinizing Hormone) ELISA Kit	Sandwich	Mouse	0.31~20ng/mL	0.19ng/mL	19,800	9,900
LH	E-EL-R0026	Rat LH(Luteinizing Hormone) ELISA Kit	Sandwich	Rat	1.56~100mIU/mL	0.94mIU/mL	19,800	9,900
LL-37	E-EL-H2438	Human LL-37(Antibacterial Protein LL-37) ELISA Kit	Sandwich	Human	1.56~100ng/mL	0.94ng/mL	24,360	12,180
MCP-1	E-EL-H6005	Human MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	Sandwich	Human	62.5~4000pg/mL	37.5pg/mL	24,360	12,180
MDA	E-EL-0060	MDA(Malondialdehyde) ELISA Kit	Competitive	Universal	31.25~2000ng/mL	18.75ng/mL	19,800	9,900
MR-ProADM	E-EL-H5441	Human MR-ProADM(Mid-regional pro-adrenomedullin) ELISA Kit	Sandwich	Human	15.63~1000pg/mL	9.38pg/mL	24,360	12,180
MT	E-EL-H2016	Human MT(Melatonin) ELISA Kit	Competitive	Human	15.63~1000pg/mL	9.38pg/mL	19,800	9,900
NA/NE	E-EL-0047	NA/NE(Noradrenaline/Norepinephrine) ELISA Kit	Competitive	Universal	0.31~20ng/mL	0.19ng/mL	19,800	9,900
NGAL	E-EL-H0096	Human NGAL(Neutrophil Gelatinase Associated Lipocalin) ELISA Kit	Sandwich	Human	0.16~10ng/mL	0.1ng/mL	24,360	12,180
NT-proBNP	E-EL-H0902	Human NT-proBNP(N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit	Sandwich	Human	0.63~40ng/mL	0.38ng/mL	24,360	12,180
NT-proBNP	E-EL-M0834	Mouse NT-proBNP(N-terminal pro-Brain Natriuretic Peptide) ELISA Kit	Sandwich	Mouse	15.625~1000pg/mL	9.38pg/mL	24,360	12,180
OC/BGP	E-EL-M0864	Mouse OC/BGP(Osteocalcin) ELISA Kit	Sandwich	Mouse	0.391~25ng/mL	0.23ng/mL	19,800	9,900
OC/BGP	E-EL-R0243	Rat OC/BGP(Osteocalcin) ELISA Kit	Sandwich	Rat	0.78~50ng/mL	0.47ng/mL	19,800	9,900
OPG	E-EL-H1341	Human OPG(Osteoprotegerin) ELISA Kit	Sandwich	Human	0.16~10ng/mL	0.1ng/mL	19,800	9,900
OT	E-EL-0029	OT(Oxytocin) ELISA Kit	Competitive	Universal	15.63~1000pg/mL	9.38pg/mL	19,800	9,900
OxLDL	E-EL-H6021	Human OxLDL(Oxidized Low Density Lipoprotein) ELISA Kit	Sandwich	Human	62.5~4000pg/mL	37.5pg/mL	19,800	9,900
P1NP	E-EL-H0185	Human P1NP(Procollagen I N-Terminal Propeptide) ELISA Kit	Sandwich	Human	15.63~1000pg/mL	9.38pg/mL	24,360	12,180
P1NP	E-EL-M0233	Mouse P1NP(Procollagen I N-Terminal Propeptide) ELISA Kit	Sandwich	Mouse	93.75~6000pg/mL	56.25pg/mL	24,360	12,180
PIIINP	E-EL-H0183	Human PIIINP(N-Terminal Procollagen III Propeptide) ELISA Kit	Sandwich	Human	23.44~1500pg/mL	14.06pg/mL	19,800	9,900
PAP	E-EL-H2101	Human PAP(Plasmin-Antiplasmin Complex) ELISA Kit	Sandwich	Human	3.13~200ng/mL	1.88ng/mL	24,360	12,180
PGE2	E-EL-0034	PGE2(Prostaglandin E2) ELISA Kit	Competitive	Universal	31.25~2000pg/mL	18.75pg/mL	19,800	9,900
PKM2	E-EL-H6063	Human PKM2(Pyruvate Kinase, Muscle) ELISA Kit	Sandwich	Human	15.63~1000IU/mL	9.38IU/mL	24,360	12,180
PYY	E-EL-H1237	Human PYY(Peptide YY) ELISA Kit	Sandwich	Human	31.25~2000pg/mL	18.75pg/mL	24,360	12,180
S100B	E-EL-H1297	Human S100B(S100 Calcium Binding Protein B) ELISA Kit	Sandwich	Human	31.25~2000pg/mL	18.75pg/mL	19,800	9,900
SOD1	E-EL-H1113	Human SOD1(Superoxide Dismutase 1, Soluble) ELISA Kit	Competitive	Human	62.5~4000pg/mL	37.5pg/mL	19,800	9,900
ST/5-HT	E-EL-0033	ST/5-HT(Serotonin/5-Hydroxytryptamine) ELISA Kit	Competitive	Universal	15.63~1000ng/mL	9.38ng/mL	19,800	9,900

Elabscience TOP100 ELISA kit

Target	Catalog	Product Name	Method	Reactivity	Detection Range	Sensitivity	定價	買一送一
TGF-β1	E-EL-H0110	Human TGF-β1(Transforming Growth Factor Beta 1) ELISA Kit	Sandwich	Human	31.25~2000pg/mL	18.75pg/mL	19,800	9,900
TIMP-1	E-EL-H0184	Human TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	Sandwich	Human	0.16~10ng/mL	0.1ng/mL	19,800	9,900
TL1A	E-EL-MK1946	Monkey TL1A(Tumor Necrosis Factor Related Ligand 1A) ELISA Kit	Sandwich	Monkey	0.16~10ng/mL	0.09ng/mL	27,320	13,660
TNF-α	E-EL-H0109	Human TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit	Sandwich	Human	7.81~500pg/mL	4.69pg/mL	19,800	9,900
TNF-α	E-EL-M0049	Mouse TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit	Sandwich	Mouse	31.25~2000pg/mL	18.75pg/mL	19,800	9,900
TNF-α	E-EL-R0019	Rat TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit	Sandwich	Rat	78.13~5000pg/mL	46.88pg/mL	19,800	9,900
TSH	E-EL-R0976	Rat TSH(Thyroid Stimulating Hormone) ELISA Kit	Sandwich	Rat	1.25~80ng/mL	0.75ng/mL	19,800	9,900
VEGF-A	E-EL-H0111	Human VEGF-A(Vascular Endothelial Cell Growth Factor A) ELISA Kit	Sandwich	Human	31.25~2000pg/mL	18.75pg/mL	19,800	9,900
VF	E-EL-H1763	Human VF(Visfatin) ELISA Kit	Sandwich	Human	0.31~20ng/mL	0.19ng/mL	24,360	12,180
Zonulin	E-EL-H5560	Human Zonulin ELISA Kit	Sandwich	Human	0.78~50ng/mL	0.47ng/mL	24,360	12,180
β-EP	E-EL-H0572	Human β-EP(Beta-Endorphin) ELISA Kit	Competitive	Human	15.63~1000pg/mL	9.38pg/mL	19,800	9,900